


Biotechnology Project Lab


Giovanna Gambarotta & Isabella Tarulli

The lecture is about to begin....


- Summary of the activity of the previous lesson
- Overview of the next lessons
- Welcome test
- RT-PCR for full length cloning or gene expression analysis

Identification of specific domains







3' end of the EGF-like domain + 5' end of the transmembrane domain





5' and 3' ends of the EGF-like domain




5' and 3' ends of the EGF-like domain








5' and 3' end of the EGF-like domain








cagaacggag caacctggtg aacatagcga atgggcctca ccacccaaac


851

NRG1 domains.annhyb - AnnHyb 4.946 April 14, 2012 Edit Oligo Sequence Multiple alignment Help (c) Olivier Friard 1997 Oligo | Oligo list | Sequence | Sequence list | Multi-alignment | Multi-alignment list Oligo U02322 NRG1 | beta4; Long name U02324 NRG1 I alpha2c Name U02322 NRG1 I beta4a Description AF194993 NRG1 II beta1a Length 1989 Type DNA checksum 5488 Author 04/10/2020 AF194439 NRG1 III-alpha2a Sequence header AF194438 NRG1 III-beta1a DQ176766 NRG1 III-beta3 Multiple alignment Sequence annotations Name init end comments 1 99 type I EGF-like 502 637 beta 638 691 atgtetgage geaaagaagg cagaggeaag gggaagggea agaagaagga 51 coggggated ogegggaage cogggeeege ogagggegae cogageeea<mark>g</mark> 101 cactgcctcc cagattgaaa gaaatgaaga gccaggagtc agctgcaggc 151 tocaagetag tgeteeggtg egaaaceage teegagtaet ceteaeteag 201 attcaaatgg ttcaagaatg ggaacgagct gaaccgcaaa aataaaccag 251 aaaacatcaa gatacagaag aagccaggga agtcagagct tcgaattaac 301 aaagcatccc tggctgactc tggagagtat atgtgcaaag tgatcagcaa 351 gttaggaaat gacagtgcct ctgccaacat caccattgtt gagtcaaacg 401 agttcatcac tggcatgcca gcctcgactg agacagccta tgtgtcctca 451 gagtetecca ttagaatete agttteaaca gaaggegeaa acaettette 501 atccacatca acatccacga ctgggaccag ccatctcata aagtgtgcgg 551 agaaggagaa aactttctgt gtgaatgggg gcgagtgctt cacggtgaag 601 gacctgtcaa acccgtcaag atacttgtgc aagtgcccaa atgagtttac 651 tggtgatcgt tgccaaaact acgtaatggc cagcttctac atgacttcta 701 ggaggaaaag gcaagaaaca gagaagcctc tagaaagaaa attggatcat 751 agccttgtga aagaatcgaa agcggaggaa ctctaccaga agagggtgct 801 gacaattact ggcatctgta tcgccctgct ggtggtcggc atcatgtgtg 851 togtogccta ctocaaaacc aagaagcagc ggcagaagct tcatgatcgg









Question 3	Correct	Mark 3.0 out of 3.0	P
You have to amplify ALL the following sequence, from the first to the last nucleotide; please design the FORWARD (=SENSE) primer, 10 base long and the REVERSE (=ANTISENSE) primer, 10 base long;			
5'-TGCT	ATCTG	CATGCTTAAG	GGCGCTAGGCTGACCTGAATCGTAGGCTAGCTAGCTAGCT

write only nucleotides, not symbols or numbers.

sense:

antisense:



You have to amplify ALL the following sequence, from the first to the last nucleotide; please design the FORWARD (=SENSE) primer, 10 base long and the REVERSE (=ANTISENSE) primer, 10 base long;


write only nucleotides, not symbols or numbers.


TGCTATCTGC


AGGTTATCTG antisense:

PROJECT OVERVIEW Sna B I Genebank sequence analysis pEGFP-N1 Ase I BsrG 1 (1389) ApaL I SnaB I Not 1 (1402) - Eco47 III (597) HSV TK poly A pEGFP-C3 primer design subcloning in different expression RT-PCR Nru 1 (209) vectors cloning in the vector pCRII-blunt pIRESpuro2 5.2 kb pCR®-Blunt II-TOPO* sequence analysis Xba I (625) < pA SV40 pA pCMV-Tag 4 pCMV-HA


- real time PCR
- protein quantification










The internal ribosome entry site (IRES) of the encephalomyocarditis virus (ECMV), permits the translation of two open reading frames from one messenger RNA.






Which is the difference among the recombinant proteins obtained with these three expression vectors?






PROTEIN OF INTEREST EGFP


wooclap




PROTEIN OF INTEREST

EGFP













PROTEIN OF INTEREST FLAG


PROTEIN OF INTEREST HA

PROTEIN OF INTEREST PUROMYCIN RESISTANCE


RT-PCR amplification

- 1 to clone the full length cDNA to express the protein
- 2 -to verify/quantify the expression of a gene/specific isoform


1 – where do you put primers if you want to amplify & clone the full length cDNA to express the corresponding full length protein?




1 – where do you put primers if you want to amplify & clone the full length cDNA to express the corresponding full length protein?




2 – where do you put primers if you want to verify & quantify the expression of a gene?



2 – where do you put primers if you want to verify & quantify the expression of a gene?




2 – where do you put primers if you want to verify & quantify the expression of a gene?




sometime you have genomic DNA contamination.

How can you amplify only cDNA and not contaminant genomic DNA?




Only for teaching purposes - not for reproduction or sale


1- one strategy could be to design primers on different exons separated by a big intron (≥1000 bp)




• which strategy could you use if intron is too small (<1000bp)?




• which strategy could you use if intron is too small (<1000bp)?




2 - if introns are too small (< 1000 bp) a good strategy could be to design primers across two exons



- what happens if you have a single exon?
- or if you have different isoforms sharing a single exon (and you have to put primers in the same exon to quantify all of them)?




if you design primers in a single exon,
you need to insert a specific control in your experiment,
to quantify possible genomic DNA contamination.



- What happen if you have **pseudogenes**?
- What are pseudogenes?

Pseudogenes are cDNA copies inside genomic DNA.

If you have genomic DNA contamination, you can amplify pseudogenes even if your primers are on different exons, separated by a big intron.



Only for teaching purposes - not for reproduction or sale

If you put primers on a single exon or you have pseudogenes, you have to insert a control to quantify genomic DNA contamination.

Which control?


RT

Which control can you do to verify that RNA and the reverse transcription ingredients are not contaminated by DNA (genomic, plasmidic, amplification products)?


	stock	1 sample	RT negative control	final concentration
RNA	0.1 μg/ μl	μl		1 µg/reaction
Buffer	5x	μl		1x
BSA	1 µg/µl	μl		0.1 μg/ μl
Triton	1%	μl		0.05%
dNTPs	10mM	μl		0,5mM
random primers	50 µM	μl		5 μΜ
RT	100u/µl	μl		200u/reaction
RNAsin	33u/µl	μl		33u/reaction
Water to 25µl		μl		
Tot		25 µl	25 µl	

The negative control is called....

Reverse Transcription




RT on the sample that you have to analyse



control "RT-": all the ingredients, except the reverse transcriptase (RT-) to detect genomic DNA contamination

* **IF** you find contamination, you can insert an extra control with **all ingredients except RNA** to detect contamination due to plasmidic DNA or previous amplification.



Usually this control (reagents, without RNA) is not carried out.

PCR reaction

Which control can you do to verify that the PCR reaction is not contaminated by DNA (plasmidic, amplification products)?


STOCK	RT+	RT -	PCR negative ctr	final conc
cDNA	1 µl			
5x buffer Taq	μl			-> 1x
sense primer 10 μM	μl			-> 250 nM
antisense primer 10 μM	μl			-> 250 nM
100% glycerol	μl			-> 5%
10mM dNTPs	μl			-> 100 µM
Taq polimerase 1u/ul	μl			-> 1u
H ₂ O	μl			
total	50 µl	50 µl		

PCR reaction

Which control can you do to verify that the PCR reaction is not contaminated by DNA (plasmidic, amplification products)?


STOCK	RT+	RT -	PCR negative ctr	final conc
cDNA	1 µl RT+	1 μl RT-		
5x buffer Taq	μl			-> 1x
sense primer 10 μM	μl			-> 250 nM
antisense primer 10 μM	μl			-> 250 nM
100% glycerol	μl			-> 5%
10mM dNTPs	μΙ			-> 100 µM
Taq polimerase 1u/ul	μl			-> 1u
H ₂ O	μΙ			
total	50 µl	50 µl		

- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- which control is missing?
- which result do you expect?



PCR for the gene of interest

- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- positive control





The positive control for your gene of interest is useful to be sure that the negative results are really negative results and not just PCR technical problems (that is: the primers are working well, the mix is ok).

 which control can you do to verify that the RNA is good and correctly retro-transcribed?

(that is, that the negative results are really negative results and not just technical problems...)


- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- positive control






- -You can amplify an housekeeping gene (TBP, actin, GAPDH,).
- -You expect that the housekeeping gene is expressed also by those samples that do not express your gene of interest.


- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- positive control



- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- positive control




- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- positive control




PCR for the gene of interest


- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- positive control




PCR for the gene of interest

- -#- sample
- -A- RT reaction negative control (RT-)
- -B- PCR reaction negative control
- -C- positive control







PCR for the gene of interest


PCR for the gene of interest

STOCK	Sample #	CTR A	CTR B	CTR C	MIX 4 samples	Conc. finale
cDNA	μl					
5 x buffer Taq	μl					-> 1x
Sense primer gene target 10 μM	μl					-> 200 nM
Antisense primer gene target 10 μM	μl					-> 200 nM
100% glycerol	μl					-> 5%
10mM dNTPs	μΙ					-> 200 µM
Taq polimerase 1u/ul	μΙ					-> 1u
H ₂ O	μΙ					
totale	50 µl					

PCR for the housekeeping gene (HKG)

STOCK	Sample #	CTR A	CTR B	CTR C	MIX 4 samples	Conc. finale
cDNA	μl					
5 x buffer Taq	μl					-> 1x
Sense primer HKG 10 μM	μΙ					-> 200 nM
Antisense primer HKG 10 μM	μl					-> 200 nM
100% glycerol	μl					-> 5%
10mM dNTPs	μΙ					-> 200 µM
Taq polimerase 1u/ul	μΙ					-> 1u
H ₂ O	μΙ					
totale	50 µl					

Reverse Transcription




RT on the sample that you have to analyse

control "RT-": all the ingredients, except the reverse transcriptase (RT-)

"positive control": RT on RNA extracted from a sample expressing the gene of interest


	1 "RT+"	2 "RT-"	"control +"	MIX x 3	Final concentration
RNA o H ₂ O	μl μl		_		1 µg
Buffer 5x	μl				1x
BSA 1 μg/μl	μl				0.1 μg/μl
Triton 1%	μl				0.05%
dNTPs 10mM	μl				0,5mM
Random primers 50µM	μl				5 µM
RT 100u/µl	μl	_			200u/25μl
RNAsin 33u/µl	μl				33u/25μl
water to 25µl	μl				
Total	25 µl	25 µl			
4			1	·	-



 $+49 \mu l mix PCR$

PCR


STOCK	1	2	3	4	MIX x 4	
RT o	1 µl #1	1 µl #2	1µl #3	_		
H ₂ O	_	_	_	1 µl		
5 x buffer Pfu	μl					-> 1x
Sense primer 10 μM	μl					-> 200 nM
Antisense primer 10 μM	μl					-> 200 nM
100% glicerolo	μl					-> 5%
10mM dNTPs	μΙ					-> 200 µM
Taq polimerasi 1u/ul	μl					-> 1u/50µl
H ₂ O	μΙ					
total	50 µl					



Only for teaching purposes - not for reproduction or sale


Exercizes

- (as soon as you finish, we will discuss your results)
- 1-find the length of the intron between domain EGF-like and domain alpha
- 2-find the length of the intron between domain EGF-like and domain beta




- 3-find coding exons (number and length) and introns (length) of Rattus norvegicus TATA box binding protein (Tbp) mRNA, complete cds (NM_001004198.1) to decide where to design primers.
- 4-find coding exons (number and length) and introns (length) of Rattus norvegicus glyceraldehyde-3-phosphate dehydrogenase (GAPDH) mRNA, complete cds (AF106860.2) to decide where to design primer
- 5-if you still have time, you can start to prepare primers for full length cloning, using annhyb or primer3.

 To design primers on different exons separated by a big intron (≥1000 bp) you need to know the exon-intron map and length



How do you identify exons and introns in the genomic DNA?



calculates the statistical significance. Learn more We have added a new function to Primer-BLAST that helps users design primers

Tue, 29 Sep 2020 12:00:00 EST

More BLAST news...

Web BLAST



blastx


S

translated nucleotide ▶ protein

tblastn

protein ▶ translated nucleotide





Enter organism common name, scientific name, or tax id

Search

Human

Mouse

Rat

Microbes